

# **Actimove® Manus**

**Eco Wrist Brace**



- **Stabilises the wrist**
- **Anatomical brace design**
- **Low profile**

# Actimove® Manus

## Eco Wrist Brace



Provides stabilisation of the wrist without restriction of hand function.

### Indications

- Inflammatory conditions, soft tissue injuries, and overuse injuries.
- Prevent further trauma caused by repetitive activity.

### Product Benefits and Design Features

- Distal brace edge comes under the palmar crease
  - Allows free flexion of the MCP joints
- Brace is contoured at the ulnar aspect of the hand
  - No constriction of the MCP's .
- High coverage of dorsal aspect of the hand
  - Better support to restrict wrist extension. Prevents swelling at the dorsum of hand
- Elastic material
  - Easy fit around the wrist and the arm.
- Pre-shaped and removable metal stay
  - Neutral position of wrist (moderate extension). Metal stay can be moulded for a better fit. Removal of the metal splint allows the brace to be washed.

### Application



### Ordering Information

| Size |   |         | Order Number |          | Sales Unit |
|------|---|---------|--------------|----------|------------|
|      |   |         | Left         | Right    |            |
| S    | < 16 cm   | < 18 cm | 72819-43     | 72819-47 | 1 x 1      |
| M    | 16-19 cm  | 18 cm   | 72819-44     | 72819-48 | 1 x 1      |
| L    | 19-21.5 cm  | 20 cm   | 72819-45     | 72819-49 | 1 x 1      |
| XL   | > 21.5 cm   | 20 cm   | 72819-46     | 72819-50 | 1 x 1      |