

**Gypsona® BP**  
Plaster of Paris Bandage



- Creamy plaster with excellent finish
- Quick setting
- Good conformability
- Low plaster loss

# Gypsona® BP

## Plaster of Paris Bandage

Gypsona® BP plaster of Paris bandage is made from a specially woven cloth, uniformly impregnated with the finest quality, fast setting plaster of Paris that produces an even cast with an excellent finish.

The bandage is spooled on a round core for superior bandage stability whilst handling. Moisture resistant packaging improves the protection of the bandages prior to use.

### • Rapid fracture immobilisation

The fast setting plaster of Paris used on Gypsona® BP bandages ensures that fractures can be rapidly immobilised, especially important in accident and emergency departments, operating theatres and busy fracture clinics.

### • Good conformability

The interlocked weave of this cloth makes the wet bandage more stable and conformable, to mould around awkward contours. The cloth also helps to minimise plaster loss.

### • Moisture resistant packaging

Each bandage is packed in a splash-proof flow wrap. The pack offers good protection to Gypsona® BP plaster of Paris bandages as well as being easy to handle. Each pack has application instructions printed on it.

### • Special cloth design

Edging technique along the bandage ensures minimal fraying and loose threads. This gives a smoother cast finish.



*Application of a Colles' slab*



*Completion of the plaster*

## Ordering Information

ORDER NO.	DESCRIPTION	PACK CONTENTS
<b>Slabs</b>		
71988-14	10cm x 20m 5 ply slab	1 x 1 slabs
71988-15	15cm x 20m 5 ply slab	1 x 1 slabs
71988-16	20cm x 20m 5 ply slab	1 x 1 slabs
71988-27	Emergency Splint†	1 x 1 slabs
<b>Rolls - Hospital Packs</b>		
71985-04	5cm x 2.7m bandage	1 x 48 rolls
71985-05	7.5cm x 2.7m bandage	1 x 48 rolls
71985-00	10cm x 2.7m bandage	1 x 48 rolls
71985-02	15cm x 2.7m bandage	1 x 24 rolls
71985-03	20cm x 2.7m bandage	1 x 24 rolls
71988-26	90cm x 5.5m wide roll	1 x 1 roll
<b>Rolls - Primary Care Packs</b>		
71985-18	7.5cm x 2.7m bandage*	1 x 12 individually boxed
71985-19	10cm x 2.7m bandage*	1 x 12 individually boxed
71985-20	15cm x 2.7m bandage	1 x 12 individually boxed

\* Drug Tariff

† 2 x 5ply slabs (15cm x 76cm) and 2 x Crepe Bandages (7.5cm x 2.7m)